
Lab 3 Report

Math 391 - Introduction to Modern Physics Lab Professor Wayne Lau University of Michigan



Jinyan Miao

Fall 2022

Contents

		Pa	age
3.1	Int	oduction	3
3.2	Exp	erimental setup	3
3.3	Vis	ualizing the data	4
3.4	An	lyzing the data	4
	3.4.1	Determining the stopping voltage	4
		Determining the Planck's constant	
3.5		nmary	
3.6	Ref	erences	8
3.7	Exp	eriment Data	9
3.8	Coc	le	11

Lab 3 - Photoelectric Effect

Introduction

In 1887, Heinrich Hertz discovered that light would cause the ejection of electrons from the surface of metals, so-called the photoelectric effect. [1] Einstein, aware of Planck's theory on photon energy quantization, theorized that the photoelectric effect was a result of the electrons absorbing individual quanta of light, or a photon, with energy proportional to the frequency, and he proposed that the universal relationship between the maximal energy E_{max} of the ejected electron and the frequency f of the incident light is given by the following: [2]

$$E_{max} = hf - e\phi (3.1)$$

where f is Planck's constant and $e\phi$ is the work function for the photoelectron-emitting surface. Here we note that (3.1) is valid only when we have $hf \ge e\phi$.

Phototube is a type of light sensor that operates according to equation (3.1). When incoming photons strike the photocathode, electrons are knocked out of the surface of the cathode. The electrons are then attracted to an anode hence causing a current. If one applies variable retarding potential across the phototube, then by (3.1), one can find a stopping voltage V_s such that there is no current presented in the phototube, in which case we can write:

$$eV_s = E_{max} = hf - e\phi (3.2)$$


In Lab 3 of Physics 391, we verified that, as suggested by (3.2), the stopping voltage V_s is independent of the intensity of the incident light but dependent on the frequency of the incident light. We used the Hamamatsu R414 phototube in our experiment, whose cathode material is cesium antimonide. We found that for incident light is of wavelength 405~nm, 532~nm, and 635~nm, the stopping voltage is given by $1.339\pm0.038~V$, $0.671\pm0.009~V$, and $0.391\pm0.007~V$, respectively. We also estimated the Planck's constant using equation (3.2) with our data, and the estimated Planck's constant is $(5.8776\pm0.2838)\cdot10^{-34}~m^2~kg/s$.

Experimental setup

In this experiment, we use the Hamamatsu R414 phototube. The phototube is placed in a photoelectric effect apparatus (a box-shape equipment), and is wired to a source such that retarding voltage can be applied across the phototube. A laser pointer is then slipped in to the apparatus. When the laser is turned on, it serves as a photon source causing the electron on the cathode of the phototube to be ejected. The values for photocurrents in our setup is too low for detection with normal instruments, so instead we measure the voltage across an $10000\,\Omega$ resistor connecting the anode to ground, and that voltage of the resistor is directly proportional to the current in the phototube. As we vary the retarding voltage across the phototube, both the retarding voltage and the corresponding voltage across anode resistor are recorded. We repeat this procedure for three laser pointers emitting light in three different wavelengths, $405\,nm$, $532\,nm$, and $635\,nm$. We then repeat the experiment with an optical neutral density filter slipped across the laser tube. The density filter has transmitted fraction $10^{-0.3}$. Six sets of data are collected in total, three for lasers without the filter, and three for lasers with the filter.

Visualizing the data

The recorded voltage across the anode with corresponding voltage across the cathode is displaced in the following plots for the three lasers:



We notice the small shift between the third and the forth dark green data points on the upper left in Figure 3. The multimeter recording the cathode voltage stopped working while we performed the measurement, and we swithced to a new one. The shift in the data points is most likely due to the slight calibration difference of the multimeter. When analyzing the data for the $532\,nm$ laser, for consistency, we neglect the first three data points.

Analyzing the data

3.4.1 Determining the stopping voltage

In order to find the stopping voltage V_s for each case, we fit our data points to the following function, which gives a good mathematical model for the photocurrent according to the Lab Manual:

$$i_{photo} = \frac{a(V_s - V)^2 \theta(V_s - V) + b(V_S - V)^4 \theta(V_s - V)}{1 + c(V_s - V)^2 \theta(V_S - V)}$$
(3.3)


where i_{photo} is the photocurrent directly proportional to V_{anode} , V is $V_{cathode}$, and θ is the unit step function defined by the following:

$$\theta: \mathbb{R} \to \{0,1\}$$
 $x \mapsto \begin{cases} 0 & x < 0 \\ 1 & x \ge 0 \end{cases}$

Note that the ratio V_{anode}/i_{photo} is resolved in the parameters a and b in (3.3), hence we can fit (3.3) by replacing i_{photo} as V_{anode} without affecting the best fit value for V_s . We fit the data points using the Python scipy function $optimize.curve_fit$.

First, we perform the fitting procedure for the reduced-intensity light data. For the incident light of wavelength $635 \, nm$, the best fit values and the standard deviations of the parameters are given by the followings:


	a	b	c	V_s
best fit value	257.1805	-899.2726	-12.2266	0.3887
standard deviation	21.2474	97.0962	3.7728	0.0069



Here the stopping voltage for reduced-intensity light of $\lambda = 635\,nm$ is $0.3887 \pm 0.0069\,V$.

For the reduced-intensity light of wavelength $532\,nm$, the fitted parameters are given by:


	a	b	c	V_s
best fit value	14.3785	-20.6220	-0.0516	0.6786
standard deviation	2.6529	4.9268	3.5398	0.0139



Here the stopping voltage for reduced-intensity light of $\lambda = 532\,nm$ is $0.6786 \pm 0.0139\,V$.

For the reduced-intensity light of wavelength 405 nm, the fitted parameters are given by:


	a	b	c	V_s
best fit value	46.9950	-12.6291	0.9009	1.3636
standard deviation	12.4997	9.0121	0.9864	0.0401



Here the stopping voltage for reduced-intensity light of $\lambda = 405\,nm$ is $1.3636 \pm 0.0401\,V$.

Now we perform the fitting procedure for the full-intensity light data. For the incident light of wavelength $635 \, nm$, the best fit values and the standard deviations of the parameters are given by the followings:


	a	b	c	V_s
best fit value	569.1897	-2151.0230	-14.3998	0.3912
standard deviation	43.8081	149.7486	2.4877	0.0065



Here we see that the stopping voltage for full-intensity light of $\lambda=635\,nm$ is estimated to be $0.3912\pm0.0065\,V$, which falls within one standard deviation from the estimation of that of the reduced-intensity light of the same wavelength.

For the full-intensity light of wavelength $532\,nm$, the fitted parameters are given by:


	a	b	c	V_s
best fit value	54.7024	-24.8958	2.5582	0.6714
standard deviation	15.2971	189.9723	12.172	0.0094



Here we see that the stopping voltage for full-intensity light of $\lambda=532\,nm$ is estimated to be $0.6714\pm0.0094\,V$, which also falls within one standard deviation from that of the half-intensity light of the same wavelength.

For the reduced-intensity light of wavelength $405 \, nm$, the fitted parameters are given by:

	a	b	c	V_s
best fit value	124.987	-22.9975	1.4588	1.3394
standard deviation	35.9933	26.4706	1.2834	0.0379



Here we see that the stopping voltage for full-intensity light of $\lambda = 405\,nm$ is estimated to be $1.3394 \pm 0.0379\,V$, which also falls within one standard deviation from that of the half-intensity light of the same wavelength.

In all wavelengths, we see that the difference between the estimated stopping voltage with or without the density filter is small (less than one standard deviation), which leads to the conclusion that the stopping voltage is independent on the intensity of the incident light. While on the other hand, we see that the stopping voltage increases as wavelength of the incident light decreases. Since we have $f = 1/\lambda$, this confirms that stopping voltage increases as frequency of incident light increases. This result is well predicted by equation (3.2), and we will verify the linear relationship between the two in the next section.

3.4.2 Determining the Planck's constant

In this section, we fit out data using equation (3.2) to find the Planck's constant:

$$eV_s = hf - e\phi (3.2)$$


where V_s is obtained from the last section. Here we have two parameters to fit, the Planck's constant h and the work function $e\phi$ of the cathode material. We have two sets of data, the one without the density filter, and the one with the density filter. Again, we use the Python scipy function *optimize.curve_fit* to get the results.

We start by fitting the set with the density filter. The best-fit parameters that we obtain and their standard deviation are given by the followings:

	h	ϕ
best fit value	$5.7191 \cdot 10^{-34}$	1.3123
standard deviation	$2.8751 \cdot 10^{-35}$	0.1081

For the one without the density filter, the parameters are given by the followings:

	h	ϕ
best fit value	$5.8776 \cdot 10^{-34}$	1.3613
standard deviation	$2.8384 \cdot 10^{-35}$	0.1067



The theoretical value for h is given by:

$$h = 6.6261 \cdot 10^{-34} \, m^2 \, kg/s$$

which differs from each of the experimental values that we estimated using the two data sets by more than 3 standard deviations:

$$\begin{cases} h = (5.7197 \pm 0.2875) \cdot 10^{-34} \, m^2 \, kg/s & \text{reduced-intensity} \\ h = (5.8776 \pm 0.2838) \cdot 10^{-34} \, m^2 \, kg/s & \text{full-intensity} \end{cases}$$

hence we are not able to statistically verify the theoretical value of the Planck's constant. This result is most likely due to the stray light that impinges on the phototube anode and cathode. Even though efforts have been made, such as covering the box on the gray foam pad when taking data, to avoid the effect of stray light, it is experimentally difficult to completely block the stray light because the box is not airtight. To get a more precise measurement of the Planck's constant, one might need to design the experiment in another way to minimize the effect of the stray light. Nevertheless, the results from our data, as shown in Fig. 10, do justify a linear relationship between the frequency of the incident light and the stopping voltage V_s , and this result is well predicted by equation (3.2).

Lastly, according to Toshimichi Sakata, the work function for the Cs₃Sb photo-cathode material is around $(1.8 \pm 0.1) \, eV$.^[3] In our experimental data, the work function of the photo-cathode is estimated to be:

$$\begin{cases} (1.3123 \pm 0.1081) \, eV & \text{reduced-intensity} \\ (1.3613 \pm 0.1067) \, eV & \text{full-intensity} \end{cases}$$

hence we see that the theoretical value of the work function differs from our experimental value by more than 3 standard deviation, in which case we are not able to conclude that the experimental result agrees with the theoretical value. And again, this is most like due to the effect of stray light that impinges on the phototube anode and cathode.

Summary

In this lab, intending to verify the universal relationship in the photoelectric effect given by equation (3.1) and (3.2), we found that the stopping voltage of the photocurrent is independent of the intensity of the incident light by comparing our experimental data obtained with and without the density filter. Even though we were not able to verify the theoretical value of the Planck's constant and the work function of the cathode using our experimental data, our data does confirm that the relationship between the stopping voltage and the frequency of the incident light is linear as predicted by equation (3.2).

References

- 1. Über den Einfluss des ultravioletten Lichtes auf die electrische Entladung, Heinrich Hertz, Annalen der Physik 267, 983–1000 (1887).
- 2. Über einen die Erzeugung und Verwandlung des Lichtes betreffenden heuristischen Gesichtspunkt, [On a Heuristic Viewpoint Concerning the Production and Transformation of Light], Albert Einstein, Annalen der Physik 17, 132–148 (1905).
- 3. Studies on the Cs3Sb Photo-Cathode, Toshimichi Sakata, Journal of the Physical Society of Japan 8, 723-730 (1953).

Experiment Data

green (532 nm)		
Vcathode	Vanode	Vanode0.3
0.05	3.65	0.68
0.1	3.35	0.642
0.15	2.9	0.6
0.2	2.8	0.518
0.25	2.25	0.44
0.3	1.7	0.357
0.35	1.22	0.272
0.4	0.8	0.189
0.45	0.476	0.113
0.5	0.22	0.06
0.55	0.1	0.025
0.6	0.024	0.01
0.65	0	0
0.7	-0.015	-0.004
0.75	-0.02	-0.006
0.8	-0.02	-0.006
0.85	-0.019	-0.002
0.9	-0.021	-0.006
0.95	-0.02	-0.007
1	-0.014	-0.007
1.05	-0.014	-0.007
1.1	-0.016	-0.012
1.15	-0.017	-0.01
1.2	-0.017	-0.008
1.25	-0.019	-0.011
1.3	-0.017	-0.01
1.35	-0.02	-0.006
1.4	-0.017	-0.005

red (635 nm)		
Vcathode	Vanode	Vanode0.3
0.05	7.23	3.32
0.1	6	2.74
0.15	4.14	1.87
0.2	2.27	1.03
0.25	1	0.445
0.3	0.338	0.145
0.35	0.062	0.03
0.4	-0.022	-0.01
0.45	-0.046	-0.02
0.5	-0.05	-0.021
0.55	-0.055	-0.023
0.6	-0.056	-0.024
0.65	-0.059	-0.025
0.7	-0.06	-0.026
0.75	-0.061	-0.027
0.8	-0.062	-0.028
0.85	-0.063	-0.029
0.9	-0.065	-0.029
0.95	-0.066	-0.03
1	-0.066	-0.032

blue (405 nm)		
Vcathode	Vanode	Vanode0.3
0.1	52.6	24.5
0.2	47.44	22.25
0.3	41.4	19.42
0.4	35.2	16.65
0.5	28.8	13.35
0.6	22.78	10.49
0.7	17.25	7.92
0.8	11.6	5.32
0.9	7.35	3.66
1	3.86	1.74
1.1	1.55	0.662
1.2	0.15	0.086
1.3	-0.41	-0.165
1.4	-0.56	-0.252
1.5	-0.63	-0.282
1.6	-0.68	-0.297
1.7	-0.7	-0.311
1.8	-0.725	-0.318
1.9	-0.74	-0.323
2	-0.758	-0.333
2.1	-0.77	-0.337
2.2	-0.78	-0.34
2.3	-0.79	-0.34
2.4	-0.8	-0.35

3.8. CODE

Code

The code for computing statistics of the data sets is attached.

```
1 UNIQNAME = "jmiu"
NAME = "Jinyan Miao"
COLLABORATORS = "Chi Han, Garen Dye"
5 import numpy as np
6 from matplotlib import pyplot as plt
7 import pandas as pd
8 from scipy import special
9 import scipy.optimize as opt
10 from scipy.optimize import curve_fit
12
13 # Read in data here
blue_data_405 = pd.read_csv('data/blue_data_405.csv',header=0)
red_data_635 = pd.read_csv('data/red_data_635.csv',header=0)
green_data_532 = pd.read_csv('data/green_data_532.csv',header=0)
18
19 print(blue_data_405)
y_data_b1 = blue_data_405['Vanode'].values
y_data_b2 = blue_data_405['Vanode0.3'].values
x_data_b1 = blue_data_405['Vcathode'].values +0.001*y_data_b1
zs x_data_b2 = blue_data_405['Vcathode'].values +0.001*y_data_b2
y_data_r1 = red_data_635['Vanode']
y_data_r2 = red_data_635['Vanode0.3']
x_data_r1 = red_data_635['Vcathode'].values +0.001*y_data_r1
31 x_data_r2 = red_data_635['Vcathode'].values +0.001*y_data_r2
y_data_g1 = green_data_532['Vanode']
y_data_g2 = green_data_532['Vanode0.3']
x_data_g1 = green_data_532['Vcathode'].values +0.001*y_data_g1
36 x_data_g2 = green_data_532['Vcathode'].values +0.001*y_data_g2
39 name = ['635', '532', '405']
40
_{
m 41} ## It is always good to first plot data you either generate or read in
42 plt.title('Figure 1 (405 nm light)', fontsize='xx-large')
43 plt.scatter(x_data_b1, y_data_b1, c='blue', label='full intensity')
44 plt.scatter(x_data_b2, y_data_b2, color=(0.5,0.5,1), label='filtered')
45 plt.grid()
46 plt.xlabel('$V_{cathode}$ (V)', fontsize='xx-large')
47 plt.ylabel('$V_{anode}$ (mV)', fontsize='xx-large')
48 plt.legend()
49 plt.show()
51
52 ## It is always good to first plot data you either generate or read in
plt.title('Figure 2 (635 nm light)', fontsize='xx-large')
plt.scatter(x_data_r1, y_data_r1, c='red', label='full intensity')
plt.scatter(x_data_r2, y_data_r2, color=(1,0.5,0.5), label='filtered')
56 plt.grid()
57 plt.xlabel('$V_{cathode}$ (V)', fontsize='xx-large')
58 plt.ylabel('$V_{anode}$ (mV)', fontsize='xx-large')
59 plt.legend()
60 plt.show()
_{\rm 63} ## It is always good to first plot data you either generate or read in
64 plt.title('Figure 3 (532 nm light)', fontsize='xx-large')
65 plt.scatter(x_data_g1, y_data_g1, c='green', label='full intensity')
66 plt.scatter(x_data_g2, y_data_g2, color=(0.5,1,0.5), label='filtered')
67 plt.grid()
plt.xlabel('$V_{cathode}$ (V)', fontsize='xx-large')
plt.ylabel('$V_{anode}$ (mV)', fontsize='xx-large')
70 plt.legend()
71 plt.show()
```

```
72
73
74
_{75} ## define the mathematical function in Eq 9 below using a, b, c, and vs as
      parameters
^{76} ## and np.heaviside(x1,x2), where we replace x1 and x2 with appropriate
_{78} ## reminder: to take a quantity X to the nth power in python, use X**n
79
   def Eq9_model(x, a, b, c, vs):
        ' A model to fit the data which is flat above the retarding potential
81
       and increases polynomially below the retarding potential
83
       x (array-like) : independent variable
84
       a (float) : model parameter
85
       b (float) : model parameter
86
87
       c (float)
                  : model parameter
       vs (float) : model parameter, the stopping voltage
88
89
     Returns:
     y (array-like) : model output values for each x
90
91
     return (a * np.heaviside(vs-x,1) * (vs-x)**3 + b * np.heaviside(vs-x,1) *
92
       (vs-x)**4)/(1 + c * np.heaviside(vs-x,1) * (vs-x)**2)
93
94
95
96 ii=[1,4,1]
97 jj=[15,19,18]
98 name = ['635', '532', '405']
99 colors = ['r', 'g', 'b']
100
_{101} ## "The true retarding potential across the phototube is obtained by adding
      the
      magnitudes of the potentials measured at the two BNC output jacks."
102 ##
x_{data_1} = [x_{data_r1}, x_{data_g1}, x_{data_b1}]
104 x_data_2 = [x_data_r2, x_data_g2, x_data_b2]
105 y_data_1 = [y_data_r1, y_data_g1, y_data_b1]
y_data_2 = [y_data_r2, y_data_g2, y_data_b2]
107 v_s2 = []
108
109
110 for k in range(3):
    ## use curve_fit starting with the data at index i and ending at index j
    i=ii[k]
112
113
    j=jj[k]
114
     x2 = x_{data_2}[k]
115
     y2 = y_data_2[k]
116
117
118
     params, params_covariance = opt.curve_fit(Eq9_model, x2[i:j], y2[i:j], p0
       = [1,1,1,0.4]
119
     perr = np.sqrt(np.diag(params_covariance)) ## np.diag returns the
120
      diagonal elements [0,0] and [1,1]
121
                                                    ## the standard deviation is
      the square root of the variance
122
     ## plot the data with the fitted function
123
     plt.scatter(x2, y2, label='data points', color=colors[k]) ## plot the
124
      data
     ## plot the fit line, also starting at index i and ending at index j
125
     plt.plot(x2[i:j], Eq9_model(x2[i:j], a=params[0], b=params[1], c=params
126
      [2], vs=params[3]), 'black', label='Fit')
     plt.plot(np.full(shape=len(x2), fill_value=params[-1].round(4)), y2, color
127
       ='gray', label='$V_s$ = %.4f'%params[-1].round(4))
     plt.grid()
     plt.xlabel('$V_{cathode}$ (V)', fontsize='xx-large')
129
     {\tt plt.ylabel('\$V_{anode})\$ (mV)', fontsize='xx-large')}
130
     plt.title('Figure '+str(k+4) + ' ('+str(name[k])+' nm light, d = 0.3)',
131
      fontsize='xx-large')
     plt.annotate('fit V$_s$ st. dev. = %.4f '%perr[-1], (0.52,0.25),
132
                   xycoords='figure fraction',
133
                   fontsize='x-large')
134
```

3.8. CODE

```
plt.legend(loc='best',fontsize='large')
     plt.show()
136
     v_s2.append(params[3])
137
     print("The stopping voltage is ", params[3].round(4), "+/-", perr[3].round
138
       (4), "V")
     print("c ", params[2].round(4), "+/-", perr[2].round(4))
139
     print("b ", params[1].round(4), "+/-", perr[1].round(4))
print("a ", params[0].round(4), "+/-", perr[0].round(4))
140
141
142
143
144 ii=[1,4,1]
145 jj=[15,19,18]
146 name = ['635', '532', '405']
147 colors = ['r', 'g', 'b']
149 ## "The true retarding potential across the phototube is obtained by adding
_{\rm 150} ## magnitudes of the potentials measured at the two BNC output jacks."
x_{data_1} = [x_{data_r1}, x_{data_g1}, x_{data_b1}]
152 x_data_2 = [x_data_r2, x_data_g2, x_data_b2]
y_data_1 = [y_data_r1, y_data_g1, y_data_b1]
y_data_2 = [y_data_r2, y_data_g2, y_data_b2]
155 \text{ v s1} = []
156
157 for k in range(3):
    ## use curve_fit starting with the data at index i and ending at index j
158
159
     i=ii[k]
160
     j=jj[k]
161
162
     x1 = x_data_1[k]
     y1 = y_data_1[k]
163
164
165
     params, params_covariance = opt.curve_fit(Eq9_model, x1[i:j], y1[i:j], p0
       = [1,1,1,0.4])
166
     perr = np.sqrt(np.diag(params_covariance)) ## np.diag returns the
167
       diagonal elements [0,0] and [1,1]
                                                        ## the standard deviation is
168
       the square root of the variance
169
     ## plot the data with the fitted function
170
171
     plt.scatter(x1, y1, label='data points', color=colors[k]) ## plot the
       data
172
     \mbox{\tt \#\#} plot the fit line, also starting at index i and ending at index j
173
     plt.plot(x1[i:j], Eq9_model(x1[i:j], a=params[0], b=params[1], c=params
       [2], vs=params[3]), 'black', label='Fit')
     plt.plot(np.full(shape=len(x1), fill_value=params[-1].round(4)), y1, color
174
       ='gray', label='$V_s$ = %.4f'%params[-1].round(4))
     plt.grid()
     plt.xlabel('$V_{cathode}$ (V)', fontsize='xx-large')
176
     plt.ylabel('$V_{anode}$ (mV)', fontsize='xx-large')
177
     plt.title('Figure '+str(k+7) + ' ('+str(name[k])+' nm light, d = 0)',
178
       fontsize='xx-large')
     plt.annotate('fit V$_s$ st. dev. = %.4f '%perr[-1], (0.52,0.25),
179
                    xycoords='figure fraction',
180
                    fontsize='x-large')
181
182
     plt.legend(loc='best',fontsize='large')
     plt.show()
183
184
     v_s1.append(params[3])
     print("The stopping voltage is ", params[3].round(4), "+/-", perr[3].round
185
       (4), "V")
     print("c ", params[2].round(4), "+/-", perr[2].round(4))
print("b ", params[1].round(4), "+/-", perr[1].round(4))
print("a ", params[0].round(4), "+/-", perr[0].round(4))
186
187
188
189
190
\ensuremath{^{191}} ## space to make a linear fit of the stopping voltage as a function of
       incident light frequency
192
193 def Eq10toFit(f, h, phi):
194
       f: frequency of light
196
       h: planck's constant
197
phi: workfunction
```

```
output:
    e*v_s = h*f - e*phi
200
201
   return h*f - phi
202
203
204 e = 1.60217663e-19
205 c=299792458
206 f = [c/(635e-9), c/(532e-9), c/(405e-9)]
207 params1, params_covariance1 = opt.curve_fit(Eq10toFit, np.array(f)/e, np.
    array(v_s1), p0 = [1, 1])
208 params2, params_covariance2 = opt.curve_fit(Eq10toFit, np.array(f)/e, np.
    array(v_s2), p0 = [1, 1])
209
perr1 = np.sqrt(np.diag(params_covariance1))
211 perr2 = np.sqrt(np.diag(params_covariance2))
212
213
214 print(params1, perr1)
215 print(params2, perr2)
216
217
plt.title('Figure 10', fontsize='xx-large')
plt.scatter(f, np.array(v_s1), c='orange', label='data (full I)')
[0])*10**34))
223 plt.grid()
plt.xlabel('frequency of light (Hz)', fontsize='xx-large')
plt.ylabel('$V_s$ (V)', fontsize='xx-large')
226 plt.legend()
227 plt.show()
```